

Driver PL2303 Prolific v1.8.19.exe

- 1) Faça o Download do Driver de Comunicação em www.cefise.com.br – Produtos – X_Conversor Serial/USB – Downloads – Driver PL2303 e FTDi – Driver PL2303 prolific DriverInstaller v1.8.19.exe



The screenshot shows the CEFISE website's 'PRODUTOS' page. On the left is a vertical menu with categories: Beat Test, Cicloergometria, Dinamometria, Dinamômetro Isocinético, Ergometria, Fococélulas, Natação, Potência, Salto Vertical, Ventilometria, X_Conversor Serial/USB, and Driver PL2303 e FTDi. The 'Driver PL2303 e FTDi' category is selected. To the right, there is an image of a USB cable and a small electronic device. Below the image is a blue 'Downloads' button. Underneath the button, two download links are listed: 'PL2303 Prolific DriverInstaller v1.8.19.exe' and 'CDM-v2.12.06-WHOL-Certified.zip'. To the right of the image, there is a small text block: 'Para utilizar o Conversor Serial/USB faça o download do Driver do Conversor Serial/USB PL2303 compatível com XP e Win7. O Driver FTDi é destinado apenas as interfaces modelos x_win10'.

- 2) Execute o aplicativo **PL2303 prolific DriverInstaller v1.8.19.exe**
- 3) Siga as instruções das caixas de diálogo.